Infratec[™] 1241 Grain Analyser



InfratecTM 1241 whole grain analyser for automatic, fast and accurate analysis of grain and flour quality parameters.

Features and benefits

- Artificial Neural Network calibrations for unique performance
- High performance scanning spectroscopy for robustness and consistency
- Truly transferable calibrations reduces running costs
- Automatic multigrain cell saves time and eliminates human errors
- Independent of sample temperature for accuracy
- Capability for high moisture samples for a wide range of applications
- Datalogger / Datalink for easy transfer of results to LIMS
- Remote control for automated analysis

Infratec™ 1241 Grain Analyser

Infratec 1241 is a dedicated whole grain analyser using Near-Infared transmittance technology to analyse a wide range of parameters (moisture, protein, oils, starch etc.) in an even wider range of commodities. The extended wavelength range of 570 - 1100 nm (standard since 2007) opens up possibilities for measurement of colour and colour-related parameters. The built-in multi-grain sample cell automatically adjusts pathlength to the corresponding application and eliminates the time, labour and risk of error otherwise related to manually adjusting pathlengths.

The analyser can be expanded with optional modules:

- · Test Weigh Module for measurement of volume weight
- Sample Transport Module for liquids, moist samples and small sample volumes
- Flour Module for measuring on flour, semolina, soymeal and other ground samples.

Global standard for grain trading

Infratec 1241 is the solution for all steps in the agricultural chain, from crop management, grain trading, and grain processing to flour production. Analyse quality parameters such as protein, moisture, starch, oil, volume weight, colour etc. with high accuracy. The Infratec 1241 Grain Analyser can be used as a stand-alone or as

The Infratec 1241 Grain Analyser can be used as a stand-alone or as a networked unit via MOSAIC software. Transfer of data between instrument and LIMS systems as well as remote control of the unit is facilitated by Datalogger and Datalink software packages.

Applications and calibration techniques

Infratec analysers have become the most widely used system in grain quality control due to its performance, i.e. a combination of instrumentation and calibrations. FOSS' introduction of ANN (Artificial Neural Network) calibration technology has revolutionized grain testing by delivering simplicity of use in combination with unsurpassed accuracy. FOSS global ANN calibration models are today used by all major grain producing countries. The largest FOSS ANN model contains over 50.000 reference data sets, giving stable, transferable calibrations and accurate analytical results. Applications (ANN and PLS) are available covering a wide range of commodities and parameters. Contact FOSS to find out what the Infratec can do for you.



Optional Modules

Flour Module Test Weight Module Sample Transport Module

Support Software

InfratecTM File Tool, 1241

WinISITM 4, Calibration Development Software

ODIN, Application Model Maker

Infratec Scan Predictor

Infratec DataLogger (included with instrument)

FOSS DataLink

Analysis time:

MOSAIC internet network software

For support and administration of Infratec systems operated in networks, contact FOSS Analytical for further information.

Operation Data Operating Programs

Software: Menu driven

Regression programs: ANN (Artificial Neural Network);

PLS (Partial Least Squares)

No. of sub-samples 1 - 20

Sample handling and result presentation

40 seconds for 7 sub-samples 50 seconds for 10 sub-samples

Path length: Variable cell automatically controlled

from 6 - 33 mm

Result report: Presented on the display as default.

Can be sent to PC/LIMS and the

printer port

Outlier function: Warnings and options for the presen-

tation of the result

Technical data

Voltage: 220-240V 50-60Hz or 110-120V

50-60Hz

Rated current: 1.0A (110-120V) / 0.5A (220-240V)

Dimensions: $W \times D \times H - 500 \times 570 \times 400 \text{ mm}$

Weight 30 kg

Monochromator: Scanning
Wavelength range: 570 - 1100 nm

Optical bandwidth: 7 nm

Number of data

points/scan: 265

Mode: Transmittance

Light source: Tungsten halogen lamp

Detector: Silicon

Storage Media: Flash disk, USB memory stick

Display: $640 \times 480 \text{ TFT LCD}$

Interface

Printer: 25 pin parallel port Modem: 9 pin serial port External PC: 9 pin serial port

LAN: RJ45 Keyboard/Barcode: PS/2 USB Ports: 2 pcs

Remote I/O: 15-pin High Density DSUB

Diagnostics: Self tests for internal communication,

monochromator and detector (offset,

gain and noise)

System protection: Dust and humidity protected

PATENTED METHOD - US PATENTS; US 4,944,589 AND EUROPEAN PATENTS; EP 0 320 477 B1, 8704886-4.

FOSS

FOSS Analytical Slangerupgade 69 DK-3400 Hilleroed Denmark

Tel.: +45 7010 3370 Fax: +45 7010 3371

info@foss.dk www.foss.dk

P/N 99991677, Issue 3GB, January 2011